

WHEATON

Buy 2 Cases of WHEATON CryoELITE® Cryogenic Vials and Receive 1 FREE Case*

CryoELITE Cryogenic Vials

- > **2D Data Matrix Bar Code Bottom Insert** provides a unique identifier for traceability.
- > **Unrivalled Cap Seal** exceeds DOT and IATA regulations ensuring ultimate protection of samples during transportation and demanding freeze-thaw handling.
- > **Loctagon™ Vial Skirt** provides stability in the freestanding position; locks into the CryoELITE Benchmate Rack to provide easy open and close with one hand.
- > **E-Z Ex-Traction® Conical Well Design** allows precise sample retrieval for both automated and manual systems.
- > **7 Cap Colors** for convenient sample identification.
- > **Lot Certified RNase/DNase and Endotoxin Free**, sterile and non-pyrogenic.



*Integrity Means Everything
For Your Precious Samples*

LIMITED TIME OFFER Redemption Coupon

Offer Valid thru December 31, 2011

Complete the redemption coupon below and send or mail/fax/email it along with your proof of purchase to:

WHEATON Science Products
Attention: Fulfillment Coordinator
1501 North 10th Street | Millville, NJ 08332-2038 USA
Email: joel.torres@wheaton.com
800.225.1437 (T) | 856.825.9035 (F)

Please send my **FREE Case of CryoELITE Cryogenic Vials** Cat. No. _____.

Please enclose or fax a copy of your proof of purchase: invoice, packing slip, or purchase order

Name: _____

Company: _____

Street: _____

City: _____ State: _____

Zip/Postal: _____ Country: _____

Tel: _____

E-Mail: _____

This offer is void where prohibited by law or company policy. All information must be completed to be a valid entry. Offer valid in North America only thru 12/31/2011. Limited to two redemptions per customer. *Free case must be of equal or lesser value. Not for resale.

CryoELITE Vials, Freestanding, External Thread

Wheaton Cat. No.	Size mL	Color	Writing Patch	Sterile	Dimensions mm	Qty
W985862	1.2	Natural	Yes	Yes	12 x 37	500
W985852	2	Natural	Yes	No	12 x 49	1000
W985853	2	Natural	No	No	12 x 49	1000
W985872	2	Natural	Yes	Yes	12 x 49	500
W985863	2	White	<input type="checkbox"/> Yes	Yes	12 x 49	500
W985864	2	Red	<input checked="" type="checkbox"/> Yes	Yes	12 x 49	500
W985865	2	Pink	<input checked="" type="checkbox"/> Yes	Yes	12 x 49	500
W985866	2	Yellow	<input checked="" type="checkbox"/> Yes	Yes	12 x 49	500
W985867	2	Green	<input checked="" type="checkbox"/> Yes	Yes	12 x 49	500
W985868	2	Blue	<input checked="" type="checkbox"/> Yes	Yes	12 x 49	500
W985869	3	Natural	Yes	Yes	12 x 63	500
W985870	4	Natural	Yes	Yes	12 x 78	500
W985871	5	Natural	Yes	Yes	12 x 93	500

CryoELITE Vials, Freestanding, Internal Thread

Wheaton Cat. No.	Size mL	Color	Writing Patch	Sterile	Dimensions mm	Qty
W985915	1.2	Natural	Yes	Yes	12 x 40	500
W985902	2	Natural	Yes	No	12 x 50	1000
W985903	2	Natural	No	No	12 x 50	1000
W985922	2	Natural	Yes	Yes	12 x 50	500
W985916	2	White	<input type="checkbox"/> Yes	Yes	12 x 50	500
W985917	2	Red	<input checked="" type="checkbox"/> Yes	Yes	12 x 50	500
W985918	2	Pink	<input checked="" type="checkbox"/> Yes	Yes	12 x 50	500
W985919	2	Yellow	<input checked="" type="checkbox"/> Yes	Yes	12 x 50	500
W985920	2	Green	<input checked="" type="checkbox"/> Yes	Yes	12 x 50	500
W985921	2	Blue	<input checked="" type="checkbox"/> Yes	Yes	12 x 50	500
W985923	3	Natural	Yes	Yes	12 x 63	500
W985924	4	Natural	Yes	Yes	12 x 77	500
W985925	5	Natural	Yes	Yes	12 x 91	500

To order contact your local laboratory products distributor.



WHEATON Science Products

a division of WHEATON Industries Inc.

1501 North 10th Street, Millville, NJ 08332-2038
800.225.1437 (U.S. & Canada) • 856.825.1100 • 856.825.1368 (F)

www.wheatonsci.com

Lit. No. 8715 07/2011