

WHEATON

Buy 2 Cases of **WHEATON Liquid Scintillation Vials** and Receive **1 FREE Case***

Liquid Scintillation Vials

- > 20mL size
- > Made from Wheaton 180 low potassium borosilicate glass
- > Background counts are consistent and low; ultraviolet transmission is high
- > Available with cap attached

Cat. No.	Cap Size	Cap Material	Liner Material	Qty/Case
----------	----------	--------------	----------------	----------

Vials with Screw Caps Attached

986540	22-400	Polypropylene	Foamed Polyethylene	500
986541	22-400	Polypropylene	Metal Foil	500
986542	22-400	Urea	Metal Foil	500
986546	22-400	Urea	Poly Seal Cone	500
986548	22-400	Urea	Polyethylene Disc	500
986560	24-400	Polypropylene	Foamed Polyethylene	500
986561	24-400	Polypropylene	Metal Foil	500
986562	24-400	Urea	Metal Foil	500
986568	24-400	Urea	Polyethylene Disc	500



LIMITED TIME OFFER Redemption Coupon

Offer Valid thru December 31, 2011

Complete the redemption coupon below and send or mail/fax/email it along with your proof of purchase to:

WHEATON Science Products
Attention: Fulfillment Coordinator
1501 North 10th Street | Millville, NJ 08332-2038 USA
Email: joel.torres@wheaton.com
800.225.1437 (T) | 856.825.9035 (F)

Please send my **FREE Case of Liquid Scintillation Vials** Cat. No. _____

Please enclose or fax a copy of your proof of purchase: invoice, packing slip, or purchase order.

Name: _____

Company: _____

Street: _____

City: _____ State: _____

Zip/Postal: _____ Country: _____

Tel: _____

E-Mail: _____

This offer is void where prohibited by law or company policy. All information must be completed to be a valid entry. Offer valid in North America only thru 12/31/2011. Offer limited to two redemptions per customer. *Free case must be of equal or lesser value. Not for resale.



 **WHEATON** Science Products
a division of **WHEATON** Industries Inc.

1501 North 10th Street, Millville, NJ 08332-2038
800.225.1437 (U.S. & Canada) • 856.825.1100 • 856.825.1368 (F)
www.wheatonsci.com