

WHEATON | SingleScan™ Bar Code Reader

Buy the SingleScan Bar Code Reader and Receive 1 FREE Case of CryoELITE™ Bar Code Inserts

WHEATON | SingleScan Bar Code Reader

The WHEATON SingleScan Bar Code Reader makes sample traceability easy. The SingleScan Bar Code Reader is a sophisticated 2D Data Matrix and 1D Linear Bar Code Reader capable of reading bar codes from a wide range of manufacturers and label printers.

- > Plug and play design enables easy set-up with no software installation required.
- > USB interface allows easy connection to computer.
- > Reads 1D linear bar codes and 2D data matrix bar codes on any vial or ampule.
- > Inputs decoded bar code ID into any software application where cursor is placed.
- > One year warranty.

| Cat. No. | Description | Qty |
|----------|----------------------------|-----|
| W986000 | SingleScan Bar Code Reader | 1 |

WHEATON | CryoELITE 2D Data Matrix Bar Code Bottom Insert

- > 2D Data Matrix Bar Code Inserts allow for immediate bar coding of your samples. The insert can also be applied to the vial at a future date, which eliminates jeopardizing the integrity of your sample by transferring it to another vial.
- > Insert press fits and locks into place in bottom of vial.

| Cat. No. | Size (mL) | Color | Writing Patch | Sterile | Dimensions (mL) | Qty |
|---|-----------|-------|---------------|---------|-----------------|-----|
| CryoELITE 2D Data Matrix Bar Code Bottom Inserts | | | | | | |
| W985880 | | | | Yes | | 500 |
| W985881 | | | | No | | 500 |

WHEATON Science Products
 a division of WHEATON Industries Inc.
 1501 North 10th Street, Millville, NJ 08332-2038
 800.225.1437 (U.S. & Canada) • 856.825.1100 • 856.825.1368 (F)
www.wheatonsci.com



LIMITED TIME OFFER* Redemption Coupon

Offer Valid thru December 31, 2011

Complete the redemption coupon below and send or mail/fax/email it along with your proof of purchase to:

WHEATON Science Products
 Attention: Fulfillment Coordinator
 1501 North 10th Street | Millville, NJ 08332-2038 USA
 Email: joel.torres@wheaton.com
 800.225.1437 (T) | 856.825.9035 (F)

Please send my FREE Case of CryoELITE Bar Code Inserts Cat. No. _____.

Please enclose or fax a copy of your proof of purchase: invoice, packing slip, or purchase order

Name: _____

Company: _____

Street: _____

City: _____ State: _____

Zip/Postal: _____ Country: _____

Tel: _____

E-Mail: _____

*This offer is void where prohibited by law or company policy. All information must be completed to be a valid entry. Offer valid thru 12/31/2011. Not for resale.