

MSP BRAND PIPETTE CAROUSEL STAND



- **Accepts most popular pipettes**
- **Stable, sturdy base**
- **Six pipette capacity**
- **Rotates 360°**

MSP Brand's pipette carousel stand is designed to accommodate a set of six single-channel pipettes and keep them within your reach, ready for use. The freely rotating circular rack section has been painstakingly engineered to fit and secure not only the most popular pipette types, but over twenty different pipette brands and the vast majority of all pipette brands on the market today. The list includes, but is not limited to: Gilson Pipetman® (all models, old and new), Rainin® (Classic and LTS), Eppendorf Research® Plus and various models of Labnet™ and HTL™ variable volume pipettes.

With a convenient handle knob at the top, the entire carousel stand can be easily and safely carried or repositioned, even with a full load of six pipettes. The wide heavyweight base contains a hidden steel stability plate, providing a solid foundation for the carousel and eliminates any possibility of tipping. Compact (only 6.25 inches in diameter), stable and versatile, this carousel stand is a welcome addition to any research laboratory.

MedSupply Partners

3715 Atlanta Industrial Pkwy, Ste B
Atlanta, Georgia 30331
Phone: 1-877-633-7779
www.medsupplypartners.com

Ordering Information

<u>Cat. No</u>	<u>Description</u>	<u>Price</u>
62-1015-1	Rotating Carousel Pipette Stand for up to six standard single-channel pipettes. Dimensions (w x h): 6.25 x 12in, (16 x 30cm)	\$89.50

Pipetman, Rainin and Eppendorf Research are registered trademarks of their respective corporate entities: Gilson Inc, Mettler-Toledo LLC and Eppendorf AG. Labnet and HTL are property of Corning Inc.