



# CULTURE & CENTRIFUGE TUBES



- **Round Bottom Culture Tubes**
- **Conical Centrifuge Tubes**
- **High Clarity Polymers**
- **Printed Graduation & Marking Spot**
- **Sterile; RNase & DNase Free**



## Sterile Culture Tubes

These disposable sterile tubes are made of high quality polypropylene (PP) that will withstand temperatures from -190° to +120°C. They are also available in crystal clear polystyrene (PS) with a more limited range of 0° to 70°C. Caps are made of polyethylene (PE) and have the popular two position design (up/ loose for aerobic work and down/ sealed for anaerobic cultures). Polystyrene tubes may be centrifuged at up to 2000xg. Polypropylene tubes withstand higher forces, up to 4000xg.

## Sterile Centrifuge Tubes

Conical 15ml and 50ml centrifuge tubes are made of premium polypropylene for clarity and centrifugal strength. Tubes and caps may be frozen at temperatures as low as -86°C. The maximum centrifugal force is 12,000xg for the 15ml tubes and 15,000xg for the 50 ml tubes\*. Caps have a flat top surface for easy marking and labeling. Cat. No. 62-1019-5 features a "skirted" bottom that allows the tube to stand upright without a tube rack.

\*Refers to liquids w/ a specific gravity of 1.0, using appropriate rotors, cushions & adapters.

## MedSupply Partners

3715 Atlanta Industrial Pkwy, Ste B  
Atlanta, Georgia 30331  
Phone: 1-877-633-7779  
[www.medsupplypartners.com](http://www.medsupplypartners.com)



62-1022-1

62-1022-2

## Ordering Information

### Culture Tubes (supplied with dual position caps attached)

Cat. No	Capacity	Dimensions	Material	Packaging Unit	Units Per Case	Price Per Case
62-1021-1	5ml*	12x75mm	PP	25 in heat seal bag	40	\$109.00
62-1021-8	5ml*	12x75mm	PS	25 in heat seal bag	40	\$109.00
62-1021-2	5ml*	12x75mm	PP	125 in heat seal bag	8	\$99.00
62-1021-3	5ml*	12x75mm	PS	125 in heat seal bag	8	\$99.00
62-1021-4	14ml*	17x100mm	PP	25 in notched bag	40	\$209.00
62-1021-5	14ml*	17x100mm	PP	1 per sterile peel	500	\$129.00
62-1021-6	14ml*	17x100mm	PS	25 in notched bag	40	\$215.00
62-1021-7	14ml*	17x100mm	PP	1 per sterile peel	500	\$129.00

\* 12x75 tubes are graduated to 4ml. 17x100 tubes are graduated to 12ml

### Centrifuge Tubes (supplied with screw caps attached)

Cat. No	Capacity	Dimensions	Material	Packaging Unit	Units Per Case	Price Per Case
62-1019-6	15ml	17x118mm	PP	50 tubes in foam rack	10	\$129.00
62-1019-7	15ml	17x118mm	PP	50 tubes in notched bag	10	\$119.00
62-1019-8	50ml	29x115mm	PP	25 tubes in foam rack	20	\$149.00
62-1019-9	50ml	29x115mm	PP	25 tubes in notched bag	20	\$139.00
62-1019-10	50ml	29x115mm	PP	25 tubes in notched bag	20	\$139.00
62-1019-5*	50ml	29x115mm	PP	25 tubes in notched bag	20	\$159.00

### Tube Racks (for all standard size tubes)

62-1022-2	Rack for 30x15ml tubes, pack of 5		\$50.00
62-1022-1	4-Way rack for 4 x 50ml, 12 x 15ml, 32 x 1.5ml, 32 x 0.5ml tubes, rainbow pack of 5		\$35.00

\* No. 62-1019-5 is a self-standing, sterile tube