

LABeler™ Handheld Laboratory Printer



- **Lightweight, handheld operation**
- **Variety of label tapes for many tube sizes**
- **Built in cutter - no need for scissors**
- **Prints a variety of barcodes**
- **“Ultra-Low” tape for liquid nitrogen storage**

Identifying and organizing increasingly large numbers of samples has become an important objective for today's busy research laboratory. The process of writing on sample tubes & containers can be messy and unreliable. It can also be difficult to find the time to access a computer and design a layout for expensive ink jet labels. The new LABeler™ Lab Printer provides the convenience of a lab marking pen with the benefits of a large desktop printer. Simply type in the desired information using the 50+ button key pad and press the print key...just press the cut lever and your label is instantly ready for application.

With over 12 available sizes, colors & material combinations, the LABeler™ is ideal for label identification of any laboratory vessel for nearly any laboratory application. Standard label tapes can be autoclaved, boiled or even frozen to temperatures as low as -50°C. Special Ultra-Low label tapes are also available for storage in -196°C environments or for liquid nitrogen applications



Available with 9V power supply
(also runs off 6 AA batteries)



Replaceable tape cassettes
(26ft / 792cm length)



MedSupply Partners

3715 Atlanta Industrial Pkwy, Ste B
Atlanta, Georgia 30331
Phone: 1-877-633-7779
www.medsupplypartners.com

Ordering Information

Cat. No.	Description	Price
62-1031-1	LABeler™ Handheld lab printer, with 9V power supply (U.S. Plug)*	\$159.95
62-1031-2	LABeler™ Handheld lab printer, without power supply	\$119.00
62-1031-3	Power supply (9V) for LABeler™ Lab Printer (U.S. Plug)*	\$41.49
62-1031-4	26' Cassette of 6mm lab tape, white with black print	\$21.00
62-1031-5	26' Cassette of 6mm lab tape, clear with black print	\$26.00
62-1031-6	26' Cassette of 12mm lab tape, white with black print	\$23.50
62-1031-7	26' Cassette of 12mm lab tape, clear with black print	\$29.75
62-1031-8	26' Cassette of 12mm lab tape, yellow with black print	\$23.50
62-1031-9	26' Cassette of 12mm lab tape, red with white print	\$23.50
62-1031-10	26' Cassette of 12mm lab tape, blue with white print	\$23.50
62-1031-11	26' Cassette of 12mm lab tape, green with white print	\$23.50
62-1031-12	26' Cassette of 24mm lab tape, white with black print	\$29.75
62-1031-13	26' Cassette of 24mm lab tape, clear with black print	\$38.00
62-1031-14	20' Cassette of 12mm ULT (ultra-low temp) lab tape, white with black print	\$26.00
62-1031-15	20' Cassette of 24mm ULT (ultra-low temp) lab tape, white with black print	\$38.00

* For EU plug, please add "-E" to the end of the item number