

**BENCHPETTE™ SERIES**



## ***BenchPette™*** ***Adjustable Pipettes***

***Designed for comfort,  
built for accuracy.***

# BenchPette™ Series

## Single Channel Pipettes

*The BenchPette™ series is the culmination of decades of liquid handling research. Every aspect of accuracy, reliability and user comfort has been carefully reviewed in designing these high precision pipettes. The result is an advanced pipette that exhibits superior accuracy & ergonomics and is comfortable to use, even during prolonged and repetitive pipetting.*

### ● Adjustable Volume & Locking Wheel

The BenchPette™ series features a continuously adjustable, digital micrometer for volume selection and a non-mandatory locking wheel (positioned at the top of the handle for simple one-handed activation/deactivation).

### ● Superior Accuracy & Reproducibility

The finely polished stainless piston and chemically resistant, maintenance-free seals provide unsurpassed accuracy. The air-jacketed handle protects the micrometer mechanism from the inaccuracies that can be caused by thermal expansion. This enhances continued pipetting precision even after prolonged use.

### ● Full Volume Range (0.1µL-10mL)

Nine Models are available for maximum precision at virtually every conceivable volume. However, only four models are needed to cover the standard range from 1µL to 5mL.

### ● Outstanding User Comfort

State of the art technologies and materials have been designed into the BenchPette's construction, making it an exceptionally light weight, well balanced instrument. The ergonomically shaped tip ejector button is designed and contoured to minimize pressure on the thumb. Actuation forces and spring pressures have been significantly reduced.



**MSP** BRAND  
**MedSupply Partners**

3715 Atlanta Industrial Pkwy, Ste B  
Atlanta, Georgia 30331  
Phone: 1-877-633-7779  
www.medsupplypartners.com

### Product Specifications:

<u>Cat. No.</u>	<u>Volume Range</u>	<u>Accuracy</u>	<u>Reproducibility</u>
62-1016-1	0.2-2µL	±12.0-1.5%	≤6.00-0.70%
62-1016-2	1-10µL	±2.5-1.0%	≤1.50-0.40%
62-1016-3	2-20µL	±2.5-1.0%	≤1.50-0.30%
62-1016-4	5-50µL	±2.0-1.0%	≤1.50-0.25%
62-1016-5	10-100µL	±1.8-0.8%	≤0.50-0.15%

<u>Cat. No.</u>	<u>Volume Range</u>	<u>Accuracy</u>	<u>Reproducibility</u>
62-1016-6	20-200µL	±1.8-0.8%	≤0.50-0.15%
62-1016-7	100-1000µL	±1.5-0.8%	≤0.30-0.15%
62-1016-8	1-5mL	±1.2-0.6%	≤0.30-0.15%
62-1016-9	1-10mL	±2.5-0.6%	≤0.30-0.15%
62-1015-1	Carousel Stand for up to 6 pipettes		