



STORAGE RACKS & BOXES

- **Polypropylene, withstand -80°C**
- **6 types - for 0.2ml to 50ml**
- **Compact & space efficient**
- **Autoclavable**
- **Numbered tube wells**



Test Tube Racks

MSP Brand's benchtop tube racks are made of 100% polypropylene and are available in three popular configurations. No. 62-1022-1 is a 4-way rack that can be used to hold 4x 50ml, 12x 15ml or 32x microtubes (0.5, 1.5 or 2.0ml). They can be snap-linked together to form as large a rack as is needed. 62-1022-2 is a unique rack that holds 30x 15ml and 20x 50ml. 62-1022-3 is a reversible rack for 96x 1.5ml on one side and 96x 0.5ml on the other.


Freezer/Storage Boxes

Suitable for refrigerator and freezer storage as low as -80°C, these boxes have translucent lids and numbered cavities. 62-1022-4 is a PCR workstation rack/box for 96x 0.2ml (plates, strips or individual tubes). 62-1022-6 is a 10x10 format with an attached hinged lid. 62-1022-5 is the standard 9x9 format with a friction fit lid.

MedSupply Partners

3715 Atlanta Industrial Pkwy, Ste B
Atlanta, Georgia 30331
Phone: 1-877-633-7779
www.medsupplypartners.com

Ordering Information



<u>Cat. No.</u>	<u>Description</u>	<u>Dimensions</u>	<u>Color</u>	<u>Price/pk</u>
62-1022-1	4-way rack for 4x50ml, 12x15ml, 32x0.5/1.5ml	6.9x3.7x2.0in	Rainbow Pk	\$35.00/pk5
62-1022-2	Combo-rack for 30x15ml and 20x50ml	8.1x6.7x2.4in	Orange	\$50.00/pk5
62-1022-3	Reversible rack for 96x1.5ml or 96x0.5ml	9.9x4.4x0.7in	Orange	\$42.00/pk5
62-1022-4	Box rack with lid for 96x0.2ml (tubes, strips or plate)	5.0x3.3x1.2in	Orange	\$35.00/pk5
62-1022-6	Freezer box, 100x2ml, with hinged lid	5.6x5.8x2.2in	Rainbow Pk	\$35.00/pk5
62-1022-5	Freezer box, 81x2ml, with friction fit lid	5.1x5.1x2.0in	Rainbow Pk	\$35.00/pk5