




SOLUTION RESERVOIRS

For Multi-channel Pipetting



- **Polystyrene**
- **Graduated***
- **Sharp “V” bottoms**
- **Pouring spout on all 4 corners**
- **Sterile & Disposable**

MSP Brand multi-channel solution reservoirs are designed for use with all multichannel pipettes. Standard 4 channel, 8 channel, 12 channel and various special purpose pipettes are all compatible with these disposable, polystyrene reservoirs. All versions have indented “spouts” in each corner to facilitate easy, dripless pour-off of unused reagents and solutions. Bottoms are formed with a sharply tapered “V” shape to minimize waste and reagent consumption.

Whether your application calls for bulk packaging, packs of 5, or individually wrapped reservoirs, MSP Brand reservoirs are supplied sterile and environmentally packed within biodegradable bags and recycled cardboard boxes.

MedSupply Partners

3715 Atlanta Industrial Pkwy, Ste B
Atlanta, Georgia 30331
Phone: 1-877-633-7779
www.medsupplypartners.com

Ordering Information

<u>Cat. No.</u>	<u>Description</u>	<u>Capacity</u>	<u>Qty per Pack</u>	<u>Price/pk</u>
62-1012-1	Solution Reservoirs, bulk pack	25ml	100	\$49.00
62-1012-2	Solution Reservoir, sterile, individually wrapped	25ml	100	\$99.00
62-1012-3	Solution Reservoir, 5 per sterile bag	25ml	200	\$99.00
62-1012-4	Solution Reservoir, bulk pack	50ml	100	\$49.00
62-1012-5	Solution Reservoir, sterile, individually wrapped	50ml	100	\$109.00
62-1012-6	Solution Reservoir, 5 per sterile bag	50ml	200	\$109.00
62-1012-7	Solution Reservoir, divided, bulk pack	50ml	100	\$49.50
62-1012-8	Solution Reservoir, divided, sterile, individually wrapped	50ml	100	\$109.90
62-1012-9	Solution Reservoir, divided, 5 per sterile bag	50ml	200	\$109.90