

# MSP BRAND

# MICROCENTRIFUGE TUBES

Always Supplied Sterile  
RNase,  
DNase &  
Pyrogen  
Free  
At No Extra Charge

**Syringe Permeable**  
2mm circle is thinner,  
allowing through-cap injection

**Raised Locking Rim**  
seals into groove in sidewall

**Smooth Rounded Lip**  
prevents "red thumb syndrome"

**Glass-Clear Side Walls**  
maximum transparency  
BioClear™ polypropylene

**Ultra-Smooth Interior Surface**  
Nano-polished molding process

**Advanced BioClear™ Polypropylene**  
enhances pellet visualization  
and minimizes optical distortion

**Ring and Groove Closure**  
protruding ring on cap locks  
into groove in upper sidewall  
with a distinct "click"

**Frosted Marking Surface**  
on non-graduated side

**Vivid Graduations**  
at 0.1ml and every 0.25ml

**Reinforced Conical Bottom**  
withstands extreme g-forces

• **Sterile comes Standard**  
by electron beam radiation

• **Self-Standing bag**  
with gusseted bottom

• **Stop-Pops™ Supplied**  
10 Stop-Pop cap locks in every bag

Free  
Stop-Pops  
in every pack

New Self-Standing  
ZipLock Bag

**MSP Brand** microcentrifuge tubes are designed for a broad range of research applications, including enzyme digests, DNA isolation/purification, centrifugation, incubation and general sample/reagent storage. With a variety of design advances, SureSeal Tubes improve sample visualization, safety, stability and overall ease of use.

**MSP Brand** tubes are manufactured under clean room conditions, using only **BioClear™** co-polymer, an advanced polypropylene that maximizes clarity. **BioClear™** has been formulated to be crystal clear, without compromising centrifugal strength (tested up to 28,000xg) or the usable temperature range (-90°C to +121°C).

The oversized tube cap (with frosted marking surface), combined with the unique ring & groove closure, enhances the security of the liquid-tight seal, while minimizing the required locking forces for a more ergonomic opening/closing of the tube caps.

All **MSP Brand** tubes are RNase, DNase and pyrogen free.

## MedSupply Partners

3715 Atlanta Industrial Pkwy, Ste B  
Atlanta, Georgia 30331  
Phone: 1-877-633-7779  
www.medsupplypartners.com

## Ordering Information

Cat. No.	Description	Qty
62-1008-1*	1.5mL microtube, clear, sterile, self-standing bag, 10 Stop-Pops included	500
62-1008-5	1.5mL microtube, sterile, self-standing bag, 10 Stop-Pops, assorted colors	500
62-1008-13	2.0mL microtube, clear, sterile, self-standing bag, 10 Stop-Pops included	500
62-1008-14	0.5mL microtube, clear, sterile, self-standing bag,	500
62-1008-15	0.5mL microtube, sterile, self-standing bag, assorted colors	500

\* Please use last digits (-6) for blue, (-7) for red, (-8) for green, (-9) for orange, (-10) for purple, (-11) for yellow, and (-12) for opaque black. (ex: 62-1008-6 for blue)



Also available in  
these great colors!