



# Premium Cell Strainers



- **Advanced design - faster flow through**
- **Superior ventilation - prevents clogging/spills**
- **Interlocking and stackable**
- **RNase, DNase, Pyrogen and Cytotoxin free**
- **Sterile (gamma) and individually packed**
- **Designed for use with 50mL tubes**  
Adapter for 15mL tubes included

The MSP Brand Premium Cell Strainer system is a fast and easy-to-use solution for detaching, isolating and straining cells from clusters and tissues. A polypropylene frame supports a nylon mesh that is available in three color-coded porosities: 40µm(blue), 70µm(white), and 100µm(yellow). A tiered filtration process is made possible by the MSP Brand's unique interlocking stackable design.

MSP Brand's molded grip and unique geometry ensure aseptic handling and eliminate the contamination risks associated with fingers contacting the filtered suspension, mesh and/or interior surfaces. The sterile strainer can be easily removed from the special blister packaging without touching its interior surface or mesh.

The advanced frame/housing design also contributes to enhanced performance and reduced processing time, by eliminating the "air-lock" effect. Improved aeration allows air to escape while filtrate is entering the tube, preventing messy spillovers. Simply open the sterile package, attach the strainer to the top of a 50mL tube and add your cell suspension/sample. Or, attach the provided reducing adapter for use in 15ml tubes. MSP Brand produces an optimized filtrate ready for maximum flow in separation columns and cytometers.

## MedSupply Partners

3715 Atlanta Industrial Pkwy, Ste B  
Atlanta, Georgia 30331  
Phone: 1-877-633-7779  
[www.medsupplypartners.com](http://www.medsupplypartners.com)

Effective Jan 1, 2017. Subject to change without notice.



Reducing adapter  
for 15mL tubes

Cat. No.	Description	Color	Price
62-1074-1	40µm, sterile, individually wrapped, 50/pk, with 1 reducing adapter for use with 15mL tubes	Blue	\$99
62-1074-2	70µm, sterile, individually wrapped, 50/pk, with 1 reducing adapter for use with 15mL tubes	White	\$99
62-1074-3	100µm, sterile, individually wrapped, 50/pk, with 1 reducing adapter for use with 15mL tubes	Yellow	\$99
62-1065-3	Reducing Adapters for 5mL and 15mL tubes, sterile, 25/pk	- -	\$69