

MSP BRAND

Weigh Boats



- **Antistatic and biologically inert**
- **Smooth surface/rounded corners**
- **Recessed corners; precise pouring**
- **Wide, flat bottoms - resist tipping**
- **For liquids, powders or granular samples**
- **White or black for visual contrast**

MSP Brand Weigh Boats are made of antistatic polystyrene and will withstand temperatures to 90°C, as well as exposure to diluted acids, alcohols, bases and aqueous solutions. Their flexibility and smooth interior surfaces with rounded corners allow for easy decanting of powdered samples. Their flat bottoms ensure stability and prevent tipping.

Weigh Boats are available in the traditional square shape as well as in a new diamond shape for easier pouring. Originally designed as weighing containers, they have proven ideal for dispensing, storage and as disposable mixing/quick freezing trays. They are available in black or white for improved sample contrast/visibility and packaged in ultra-compact dispenser packs.

MedSupply Partners

3715 Atlanta Industrial Pkwy, Ste B
Atlanta, Georgia 30331
Phone: 1-877-633-7779
www.medsupplypartners.com

Ordering Information

Square Weigh Boats

Cat. No.	Description	Color	Price	Qty.
62-1017-1	20ml, 45x45mm	White	\$35.00	500
62-1017-2	100ml, 80x80mm	White	\$50.00	500
62-1017-3	250ml, 140x140mm	White	\$89.00	500
62-1017-4	20ml, 45x45mm	Black	\$35.00	500
62-1017-5	100ml, 80x80mm	Black	\$50.00	500
62-1017-6	250ml, 140x140mm	Black	\$89.00	500

Diamond Weigh Boats

Cat. No.	Description	Color	Price	Qty.
62-1018-1	5ml, 35x55mm	White	\$35.00	500
62-1018-2	30ml, 65x85mm	White	\$50.00	500
62-1018-3	100ml, 100x125mm	White	\$42.50	250
62-1018-4	5ml, 35x55mm	Black	\$35.00	500
62-1018-5	30ml, 65x85mm	Black	\$50.00	500
62-1018-6	100ml, 100x125mm	Black	\$42.50	250