



WESTERN BLOT BOXES



- **A Perfect Box for Every Blot**
- **Reduce Antibody Usage**
- **Hinged Lid Minimizes Evaporation**
- **For Halves, Wholes, Strips, or Minis**

No matter what size your Western blot is, there's a box that fits it perfectly. With boxes that are correctly matched to common blot sizes, you can fit more blots on shakers/rockers while keeping reagent and antibody usage to an absolute minimum. Additionally, all corners are formed at a 90° angle and the bottom surface is uniformly flat, further reducing consumption of costly reagents and antibodies. (For example, popular mini-gels require only about 4 ml.)

At a price that makes them truly disposable, all sizes are developed from crystal clear polystyrene (except B1200-7BK & B1200-15BK, which are an opaque black that protects light sensitive samples). All sizes are also equipped with an easy-to-open hinged lid that reduces evaporation.

Cat No.	Description	Applications	Minimum Volume
62-B1200-3	Micro-strip	Various	1 ml
62-B1200-5	Mini-rectangle	Half mini-protein gels	3 ml
62-B1200-7	Medium Rectangle	Mini-protein gels, Bio-Rad Protean 3 & similar	4 ml
62-B1200-9	Long Strip	Protein 3 (cut up)	2 ml
62-B1200-11	Medium Strip	Various	3 ml
62-B1200-13	Large	Cambrex PAGER Gold, Mighty Small, Novex SureLock Mini	6 ml
62-B1200-15	Square	Owl Penguin P9DS/Puffin	7 ml
62-B1200-17	Extra Large	Bio-Rad Criterion, Amersham	11 ml
62-B1200-19	XXL	Bio-Rad Criterion, Amersham	12 ml
62-B1206	Long Strip (x6)	Protein 3 (multiple)	2 ml x 6
62-B1300-8BK	Medium Square	Novex Mini Gels	4 ml

MedSupply Partners

3715 Atlanta Industrial Pkwy, Ste B
Atlanta, Georgia 30331
Phone: 1-877-633-7779
www.medsupplypartners.com

Ordering Information

Cat. No.	Dimensions	Price Pk/5
62-B1200-3	2 ⁷ / ₈ x 1 ³ / ₁₆ x ³ / ₄ in (7.3 x 3.0 x 1.9 cm)	\$16.00
62-B1200-5	2 ⁷ / ₈ x 2 x 1 ¹ / ₄ in (7.3 x 5.1 x 3.2 cm)	\$21.00
62-B1200-7	3 ¹ / ₂ x 2 ⁹ / ₁₆ x 1 in (8.9 x 6.5 x 2.5 cm)	\$28.00
62-B1200-7BK	3 ¹ / ₂ x 2 ⁹ / ₁₆ x 1 in (8.9 x 6.5 x 2.5 cm)	\$30.00
62-B1300-8BK	8.6 x 8.6 x 2.8 cm (pull off lid)	\$69.00 (10/pk)
62-B1200-9	3 ³ / ₄ x 1 ³ / ₁₆ x ⁵ / ₈ in (9.5 x 3.0 x 1.6 cm)	\$21.00
62-B1200-11	4 ³ / ₈ x 2 ¹ / ₁₆ x 1 ¹ / ₈ in (11.1 x 5.2 x 4.4 cm)	\$34.00

Cat. No.	Dimensions	Price Pk/5
62-B1200-13	4 ⁵ / ₈ x 3 ¹ / ₂ x 1 ¹ / ₈ in (11.7 x 8.9 x 2.8 cm)	\$37.00
62-B1200-13BK	4 ⁵ / ₈ x 3 ¹ / ₂ x 1 ¹ / ₈ in (11.7 x 8.9 x 2.8 cm)	\$38.00
62-B1200-15	4 ¹ / ₂ x 4 ¹ / ₂ x 1 ¹ / ₄ in (11.4 x 11.4 x 2.8 cm)	\$47.00
62-B1200-15BK	4 ¹ / ₂ x 1 ¹ / ₂ x 1 ¹ / ₄ in (11.4 x 11.4 x 2.8 cm)	\$50.00
62-B1200-17	6 x 4 x 1 ¹ / ₄ in (15.2 x 10.2 x 2.8 cm)	\$47.00
62-B1200-19	6 x 4 x 2 in (15.2 x 10.2 x 5.1 cm)	\$48.00
62-B1206	6 compartment box (11.7 x 8.9 x 2.8 cm)	\$55.00 (4/pk)