


# MSP

BRAND

## BLU-X™ AUTORADIOGRAPHY FILM



- **High Contrast**
- **Reduced Background**
- **Ideal for Chemiluminescence**
- **Publication Quality**
- **Suitable for <sup>14</sup>C, <sup>32</sup>P, <sup>33</sup>P, <sup>125</sup>I, <sup>35</sup>S**
- **Double Emulsion**

*Blu-X is a high quality, blue sensitive film that is ideal for a wide range of electrophoretic and autoradiography applications. It can be used with calcium tungstate or blue-emitting rare-earth screens.*

*Blu-X exhibits excellent sensitivity and contrast characteristics and is coated on both sides for improved consistency and overall performance. It reliably produces sharp, bright images that are easily visualized and suitable for publication.*

*Especially useful with chemiluminescence, Blu-X works well with <sup>14</sup>C, <sup>32</sup>P, <sup>33</sup>P, <sup>125</sup>I, and <sup>35</sup>S.*

*Compatible with automatic processors and manual development methods, Blu-X compares directly with films costing up to 3 times as much.*

*With a shelf life of 2 years from date of production, Blu-X is packaged in lined, light-safe dispenser boxes of 100 sheets per box.*

### MedSupply Partners

3715 Atlanta Industrial Pkwy, Ste B  
Atlanta, Georgia 30331  
Phone: 1-877-633-7779  
[www.medsupplypartners.com](http://www.medsupplypartners.com)

### Ordering Information

<u>Cat. No.</u>	<u>Description</u>	<u>Quantity</u>	<u>Price</u>
62-1027-1	Blu-X Autoradiography film, 5x7 (127x178mm)	100 Sheets	\$79.00
62-1027-2	Blu-X Autoradiography film, 8x10 (203x254mm)	100 Sheets	\$105.00